

## X Mobil



- 
- 
- 

X Twin

## X Mobil

Mobile Patient table with single side suspended, floating table top and electromagnetic locks. Motorized height adjustment for optimal patient positioning.

- Floating table top for optimum access to patient and large radiolucent exposure area
- High mobility of the table due to swivel Castors and for height adjustment by re-chargeable batteries
- Optimal in combination with the Digital Radiography System X Twin



### Mobile Patient Table

	Floating table top
1700 mm	Useful radiolucent area
2150 mm	Total table top length
650 mm	Table top width
560 mm	Table top longitudinal movement
+/- 120 mm	Table top lateral movement
Carbon	Table material
<1,1 mm AL	Attenuation equivalent
	Single button locking system for floating table top
250 kg / 550 lbs	Maximum patient weight
610 +/- 10 mm – 1010 +/- 10 mm	Table height
	Power supply

### Handling

- Double-sided locking
- Operating from both sides of the table top
- Foot pedal for lock / unlock of swivel Castors

/