

## X Mobil Q




- 
- 
- 

X Twin

## X Mobil Q

Mobile Patient table with single side suspended, floating table top and electromagnetic locks.

- Floating table top for optimum access to patient and large radiolucent exposure area
- High mobility of the table due to swivel Castors and re-chargeable battery
- Optimal in combination with the Digital Radiography System X Twin



### Mobile Patient Table

Floating table top

Useful radiolucent area

Total table top length

Table top width

Table top longitudinal movement

Table top lateral movement

Table material

Attenuation equivalent

Single button locking system for floating table top

Maximum patient weight

Table height

Power supply

weight

1700 mm

2150 mm

650 mm

560 mm

+/- 100 mm

Melamin

<1,2 mm AL

250 kg / 550 lbs

740 +/- 10 mm

Akku

150 kg

/

### Handling

Double-sided locking

Operating from both sides  
of the table top

Foot pedal for lock / unlock  
of swivel Castors